Ultimate Healthcare

4 bar side rails

Product Code: UBEB7002 (beige) Product Code: UBEB7002G (grey)



The 4-bar Side Rails are designed to fit the Cura II range and will provide sufficient height to accommodate a typical mattress /overlay combination to enhance patient safety.

Featuring a safety locking/unlocking mechanism which prevents unintentional activation, the 4-bar side rails provide an overall height of 470mm and can be folded into the desired position with ease.

Features

- Overall height of 470mm
- Safety locking/unlocking mechanism
- Easy fitting without need for tools
- Conforms to IEC60601-2-52

Specifications

Overall height (from platform):470mmOverall length:1430mm

Ultimate Healthcare Calmore Industrial Estate, Nutwood Way, Totton, Hampshire, SO40 3WW tel: 0333 321 8996, fax: 0845 051 0035, email: sales@ultimatehealthcare.co.uk



www.ultimatehealthcare.co.uk